

# QuickSeal qPCR Optical Film™

## Pressure Sensitive Seal

<b>Description</b>	An optically clear, DMSO resistant pressure sensitive seal which is suited for qPCR (96 or 384well) fluorescence, crystallation, storage. A transparent non-tacky film which adheres only when pressure is applied. It is pierceable and peelable. Good temperature and chemical resistance and withstands tough application environments. High Adhesion Strength.			
<b>Ordering</b>	IST-128-080LS IST-128-080TS	Std LabSheet™ Trial LabSheet™	100 Sheets 5 Sheets	140mm x 80mm 140mm x 80mm
<b>Application</b>	qPCR (94 or 384well) and situations where fluorescence is experienced.			
<b>Storage</b>	Store in a cool place. Avoid direct exposure to sunlight. Use seals within the expiry date shown on the label.			
<b>Properties</b>	Temperature range -80 °C to 120 °C			
<b>Sealing</b>	IST-128 Recommended sealing Equipment: KAPS 500/Sealit100/Manual Roller IST-202-102HR			